



Prospective Evaluation of the ResistancePlus™ MG:

Simultaneous Detection of *Mycoplasma genitalium* (MG) and Azithromycin Resistance Mutations

Litty Tan, PhD
Director of R&D

Investigators



Royal Women's Hospital (RWH), Melbourne, Australia

- Sepehr Tabrizi
- Jimmy Twin
- Jenny Su
- Marin Poljak
- Suzanne Garland

Melbourne Sexual Health Centre (MSHC), Australia

- Catriona Bradshaw
- Christopher Fairley

SpeedX

- Samantha Walker
- Elisa Mokany

Disclosure

SpeedX is the developer and manufacturer of the *ResistancePlus™* assay evaluated in this study

***Mycoplasma genitalium* (MG) clinical association**



General population

- § Detected in 1-3.3%
- § Exact burden of disease is unknown

Men

- § Strongly associated with non-gonococcal urethritis (NGU) (6-50%)
- § Other conditions: Proctitis, Balanitis, Chronic prostatitis, Acute epididymitis

Women

- § Strongly associated with cervicitis (5-20%)
- § Other conditions: Endometritis, Pelvic inflammatory disease, Preterm birth

3

Clinical management of MG



Empiric treatment of NGU & cervicitis


- § **Doxycycline:** poor efficacy, microbiological cure rate of MG is 30-40%
- § **Azithromycin:** higher efficacy, 1 g single dose has 85% cure rate

Decreasing overall cure rates with azithromycin due to MG resistance

- § Where 1 g azithromycin is used to treat NGU, macrolide resistance is ~30-40%

4

Macrolide resistance is an emerging global issue



Bissessor et. al., Clin. Infect. Dis. 2015;60(8):1228-36
 Chrisment et al. J. Antimicrob. Chemother. 2012;67(11):2598-601
 Touati et. al, J. Clin. Microbiol. 2014;52(5):1549-55.
 Le Roy, Emerg. Inf. Dis, 2016;22(9):1677-9.

5

Macrolide resistance in MG



Azithromycin resistance is associated with mutations in 23S rRNA gene

- § Prevent azithromycin binding and inhibiting translation
- § A2058G, A2059G, A2058T, A2058C, A2059C (*E. coli* numbering)
 - § Confer high MICs or treatment failure
- § 82% of treatment failures showed presence these mutations prior to treatment
 - § 18% showed a mutation post-treatment

Bissessor et. al., Clin. Infect. Dis. 2015;60(8):1228-36

6

MG testing



MG detection

- § In house qPCR tests and recent commercially available CE marked tests

Macrolide resistance detection (23S rRNA mutation)

- § Methods for mutation detection
 - § Sequencing – Costly and generally not convenient for routine diagnostics
 - § High resolution melt analysis (HRMA) – Separate assay to MG detection, not easy to analyse
 - § Fluorescence resonance energy transfer (FRET) – Lacking in sensitivity

7

*ResistancePlus*TM MG

Go beyond detection


CE IVD

- § Single well qPCR assay
- § Detects MG
- § Detects 5 of the most common 23S rRNA macrolide resistance mutations

8

PlexPrime™ for mutation detection

Enhanced specificity and superior multiplexing



Mutation specific **PlexPrime** AND Mutation specific **PlexZyme** detection

Patent: WO 2013/123552

9

ResistancePlus™ MG

Sensitivity and specificity



	Target	
Channel 1	MG detection (MgPa)	
Channel 2	23S rRNA mutation	A2058G
		A2059G
		A2058T
		A2058C
		A2059C
Channel 3	Internal Control	

10

ResistancePlus™ MG

Sensitivity and specificity



	Target		Sensitivity (copies/reaction)
Channel 1	MG detection (MgPa)		10
Channel 2	23S rRNA mutation	A2058G	12
		A2059G	10
		A2058T	10
		A2058C	10
		A2059C	15
Channel 3	Internal Control		n/a
	Specific against other targets		
Targets	<i>M. hominis</i> , <i>U. parvum</i> , <i>U. urealyticum</i> , <i>C. trachomatis</i> , <i>N. gonorrhoeae</i> , <i>T. vaginalis</i> , <i>S. pneumoniae</i> , <i>M. pneumoniae</i> , <i>M. gallisepticum</i> , <i>E. coli</i> , <i>E. faecalis</i> , <i>P. aeruginosa</i> , <i>K. pneumoniae</i> , HSV-1, HSV-2		

11

Prospective study: **ResistancePlus™ MG**

Sample cohort



- 1089 samples received Nov 2015-Dec 2015 at RWH (Melbourne, Australia)
 - Royal Women's Hospital, Melbourne Sexual Health Centre, and external sites
 - Symptomatic and asymptomatic males and females

	Urine/ urethral swab	Anal swab	Cervical/ vaginal swab	Sample numbers
Male	354	34	n/a	388
Female	203	2	496	701
Total	557	36	496	1089

12

Prospective study: *ResistancePlus*TM MG *M. genitalium* detection



		In house qPCR (16S rRNA)				%	95% CI
		+	-	Total			
Speedx	+	64	0	64	Sensitivity	98.5	91.7 to 99.9
	-	1	1024	1025	Specificity	100.0	99.6 to 100.0
	Total	65	1024	1089	PPV	100.0	94.4 to 100.0
					NPV	99.9	99.5 to 100.0

MG prevalence = 6.0%

- Male 10.8% (42/388)
- Female 3.3% (23/701)

13

Prospective study: *ResistancePlus*TM MG 23S rRNA mutant detection



		Sanger Sequencing				%	95% CI
		Mutant	Wild type	Total			
Speedx	Mutant	38	1	39	Sensitivity	100.0	90.8 to 100.0
	Wild type	0	25	25	Specificity	96.2	80.4 to 99.9
	Total	38	26	64*	PPV	97.4	86.5 to 99.9
					NPV	100.0	86.3 to 100.0

* Only includes MG positive samples by both methods


MG 23S rRNA mutant prevalence = 63.1%

- Male 81.0% (34/42)
- Female 30.4% (7/23)

Highly multiplexed test with excellent clinical sensitivity and specificity
Submitted for TGA Approval

14

Improved patient care with *ResistancePlus*™ MG



15

ResistancePlus™ assays in development



- § MG fluoroquinolone resistance
- § Gonorrhoea macrolide resistance

16

MG and emerging resistance



Fluoroquinolones

- § Moxifloxacin often used as second line treatment
 - 12% treatment failure in Australia
- § Resistance linked to clustered mutations in *parC*
 - High levels of resistance linked to a **G248T** mutation
 - Other *parC* mutations in the literature include **A247C**, **G259A/C/T**

17

ResistancePlus™ MG – Fluoroquinolone resistance Preliminary clinical testing



	Target	
Channel 1	MG detection (MgPa)	
Channel 2	parC mutation	G248T
Channel 3	parC mutations	G247C
		G259A
		G259T
		G259C
Channel 4	Internal Control	

MG detection: 100% sensitivity (18/18)

MG parC mutations: 91% sensitivity (10/11) & 100% specificity (7/7)



For further queries:

Litty Tan, PhD
Director of Research and Development
littyt@speedx.com.au

Please come and visit our Booth!